


Locked Nucleic Acid - LNA™



- ▶ Get the highest sensitivity and specificity in your nucleic acid recognition assays

LNA™ : Unique Technology - Unique Applications

What is a Locked Nucleic Acid?


Locked Nucleic Acids (LNA™) are a class of nucleic acid analogues in which the ribose ring is “locked” by a methylene bridge connecting the 2'-O atom with the 4'-C atom (see structure below). LNA™ nucleosides contain the six common nucleobases (T, C, G, A, U and mC) that appear in DNA and RNA and thus are able to form base-pairs according to standard Watson-Crick base pairing rules.

Oligonucleotides incorporating LNA™ have increased thermal stability and improved discriminative power with respect to their nucleic acid targets.

LNA™ can be mixed with DNA, RNA and other nucleic acid analogs using standard phosphoramidite synthesis chemistry. LNA™ oligonucleotides can easily be labeled with standard oligonucleotide tags such as DIG, fluorescent dyes, biotin, amino-linkers, etc. Thus a very high degree of freedom in the design of primers and probes exists.

LNA™	
Tm increase/monomer against DNA	2-6°C
Tm increase/monomer against RNA	3-8°C
ΔTm at single mismatch against DNA	<8°C
Compatible with standard molecular biology	Yes
Water solubility	High

The chemical structure of Locked Nucleic Acid



Why work with LNA™?

Short probes with high Tms

- ▶ Perfect for detection of short RNA and DNA targets

Increases the thermal stability of duplexes

- ▶ High-affinity binding

Increased discriminatory power

- ▶ Single base discrimination capability

Resistant to exo- and endonucleases

- ▶ High stability for *in vivo* and *in vitro* application

Increased target specificity

- ▶ Fast binding to targets

Strand invasion

- ▶ Detect “hard to access” samples

Substrate for standard enzymes, e.g. T4 PNK, T4 ligase & DNA polymerase

- ▶ Compatible with standard enzymatic processes

When should I use LNA™ probes?

In any application with a short RNA or DNA target where normal oligonucleotides do not show sufficient affinity or specificity.

Proven applications of LNA™

- ▶ microRNA research
- ▶ Real-time PCR and ProbeLibrary
- ▶ SNP genotyping
- ▶ *In situ* hybridisation
- ▶ Microarray gene expression profile
- ▶ RNAi
- ▶ Antisense oligonucleotides
- ▶ Allele specific PCR
- ▶ Cytogenetics
- ▶ DNAzymes
- ▶ Gene repair/exon skipping
- ▶ mRNA isolation
- ▶ Splice variant detection
- ▶ Comparative Genome Hybridization (CGH)


LNA™ publications:

View all LNA™ publications on www.exiqon.com/publications

The best tool for recognition and detection of RNA and DNA targets

microRNA detection by in situ hybridisation

miRCURY™ LNA Detection probes reveal new information about miRNA distribution




Specific detection of miR-206 using miRCURY™ LNA Detection probes in *in situ* hybridisation of whole mount zebrafish embryos.

Image from Wienholds et al., Science 2005 (309), 310-11.

mRNA *in situ* hybridisation

Fast and specific mRNA *in situ* hybridisation with specific LNA™ oligonucleotide in fixed cells




Improved signal and less background using a LNA™ mRNA *in situ* hybridisation probe (left picture) compared to a DNA probe (right picture).

Images from Thomsen et al., RNA 2005, (11), 1745 - 48.

Cytogenetics


Fast and specific detection of chromosomal sequences directly on intact chromosomes



Specific telomere detection using LNA™ FISH probes.
Image kindly provided by Dr. R. W. Dirks, Leiden University Medical Center, Leiden, The Netherlands.

microRNA expression profiling

miRCURY™ LNA Array capture probes give high sensitivity and specificity for small RNA targets



miRCURY™ LNA Arrays require only 1 µg of total RNA to profile miRNA. Identical miRNA profiles are produced from starting amounts of total RNA that span the range 10µg to 1µg, without miRNA enrichment. 17 different miRNAs detected in human lung total RNA (Ambion) are represented. Numbers in the top righthand of each box show the amount of total RNA used to produce each profile.

LNA™ products:

- ▶ **microRNA analysis**
miRCURY™ LNA Detection, Array, Knockdown, and Real-time PCR
- ▶ **Real-time PCR**
ProbeLibrary™ Real-Time PCR Assay System for gene expression analysis
- ▶ **In situ hybridisation**
Custom and pre-designed probes for mRNAs, small RNAs, snRNAs, and chromosomes
- ▶ **mRNA isolation**
Poly(A)+ RNA isolation using LNA™ oligo-T₂₀ capture probe
- ▶ **LNA™ Oligonucleotides**
LNA™ Oligonucleotides are available for a variety of different specialty applications and innovative products
- ▶ **Reagents**
Exiqon offers a variety of basic research reagents based on our core technologies - LNA™, AQ-Link™ and A-quencher



Patents and Trademarks

Exiqon, AQ-Link, LNA, miRCURY, ProbeLibrary and ProbeFinder are registered trademarks of Exiqon A/S, Vedbaek, Denmark. Anthraquinone photochemistry products (AQ products) are covered by patents owned by Exiqon A/S.

Locked Nucleic Acids are covered by the following patents owned by Prof. Imanishi: AU 720472, AU 742476, CA 2,283,509, USP 6,268,490 and USP 6,770,748) and by patents and patent applications owned by Exiqon A/S.

Disclaimers

Exiqon Products: Products are for research use only and not for diagnostic or therapeutic use. The products may be used only for the buyer's internal research purposes and not for commercial use. The buyer may not resell products in their original or any modified form. The purchase of products does not include or carry an implied right or license for the buyer to use such products in the provision of services to third parties and a license must be obtained directly from Exiqon A/S for such use.

Prespotted arrays: This product and its use are covered by one or more of the following patents owned by Oxford Gene Technology Limited or Oxford Gene Technology IP Limited: US 6,054,270, US 5,700,637, EP 0,373,203; Jap. 3,393, 528 and 3,386,391 and pending patents. The purchaser is licensed to practice methods and processes covered by these patents using this product for its own internal research purposes only but may not: transfer data derived from the use of this product to third parties for value; use this product in the provision of services to third parties for value; use this product to make, have made, create or contribute to the creation of stand alone expression database products for license, sale or other transfer to a third party for value; or use this product for the identification of antisense reagents or the empirical design of probes or sets of probes for using or making nucleic acid arrays.

Capture probe oligonucleotide sets: Exiqon A/S is not licensed under any patents owned by Oxford Gene Technology Limited or related companies (xOGT®) and cannot pass any such license to its customers. A license under OGT's patents may be necessary to manufacture or use oligonucleotide arrays. To enquire about a license under OGT's oligonucleotide array patents, please contact licensing@ogt.co.uk.

PCR: PCR is covered by patents owned by Roche Inc. Use of the PCR process requires a license. The use of certain probes in 5' nuclease assays is covered by European patent EP-B1-0543942 and corresponding patents. The purchase of this product does not provide a license to use the patented technology described in the above-mentioned patent. A license to practice this technology for research applications must be obtained by the end-user from the company entitled to grant licenses under this patent.

DIG: DIG is licensed from Roche Diagnostics GmbH.

Specifications in this document are subject to change without notice.